



**MedSource**  
Diagnostics

# Single and Dual Head Stethoscope



The MedSource Single Head Stethoscope is an economical choice that does not compromise on sounds quality. This unit features an aluminum chest piece, 22" Y tubing and chrome-plated brass binaurals.



The MedSource Dual Head Stethoscope, which includes the same great features as our reliable Single Head Stethoscope, but with the dual head chest piece on a 180 degree swivel. Comes complete with spare diaphragm and eartips. Available in various colors and is not made with natural rubber latex.

Rev. 1

[www.gomedsource.com](http://www.gomedsource.com)

Item #	Description	UOM
MS-70021	Single Head, Black	EA, 50 EA/CS
MS-70022	Single Head, Blue	EA, 50 EA/CS
MS-70023	Single Head, Red	EA, 50 EA/CS
MS-70025	Single Head, Yellow	EA, 50 EA/CS
MS-70031	Dual Head, Black	EA, 50 EA/CS
MS-70032	Dual Head, Blue	EA, 50 EA/CS
MS-70033	Dual Head, Red	EA, 50 EA/CS
MS-70035	Dual Head, Yellow	EA, 50 EA/CS