



MedSource
Diagnostics

Blood Glucose Meter



Glucose Meter Specifications:

- Dimensions: 94.5mm x 56mm x 27.5mm
- Weight: 56g
- Power Source: CR2032 3V Lithium coin battery x 1
- Battery Life: 1000 results
- Display: Large LCD
- External Output: RS232 PC interface
- Memory: 360 Results with time and date
- Auto electrode inserting detection
- Auto sample loading detection
- Auto reaction time count-down
- Sleep Mode: Power consumption less than 20uA
- Auto turn-off after one minute without action
- Temperature warning
- Operation Condition: 14 degree C - 40 degree C (57.2 degree F - 104 degree F)

Test Strip Specifications:

- Strip Size: 35mm x 7mm x 0.35mm
- Required Blood Sample: 1u1
- Test Result Time: 6 Seconds
- Hematocrit: 30%-55%
- Sample Type: Fresh capillary whole blood
- Test Range: 20-600mg/dL or 1.1-33.3mmol/L
- Operation Condition: 14 degree C - 40 degree C (57.2 degree F - 104 degree F)
- Validity: Can be used within 90 days after opening or within 18 months after the production date
- Unit of Packaging: 25 strips x 2 vials/box

Rev. 1

www.gomedsource.com

Item #	Description	UOM
MS-76010	Blood Glucose Meter	EA
MS-76012	Blood Glucose Test Strips	25 strips/vials, 2 vials/BX
MS-76014	Blood Glucose Lancet	100 EA/BX
MS-76016	Blood Glucose Test Solution	EA