

**SP02**  
Blood Oxygen Saturation level



**Diplays Battery Life**  
Battery low indicator

**Pulse Rate**  
Monitors heart beats per minute

**Pulse Bar & Wave Form**  
Monitors heart beats per minute

The ChoiceMMed™ Hand Held Pulse Oximeter measures the SpO2, Pulse Rate and Pulse Bar, and displays results accurately on the LED screen. The Pulse Oximeter is lightweight and portable.

Item #	Description	UOM
MS-74012	Hand Held Pulse Oximeter	EA

# Features and Benefits

- Portable, light weight Hand Held Pulse Oximeter
- Fits perfectly in your hand and weighs 165gm
- High resolution 2.8" TFT screen displays SpO2, PR wave form and pulse bar
- Adjustable audible and visual alarms
- Battery-low indicator
- (3) "AA" batteries
- (1) Adult finger probe included
- Adjustable backlight and voice
- Automatic power-off function
- 127 ID setup, 72-hour data storage and for review
- Bluetooth, charger stand and AC/DC adapter are optional
- User manual in English, Spanish, Arabic, and French
- Warranty period of 3 years

# Specifications

## Patient Range

- Adult

## SpO2

- Measurement range: 70%-100%
- Resolution: 1%
- Measurement accuracy: 70% - 100%: +3 Digit, 0%-69% no definition

## Display

- Type: 2.8" color TFT
- Resolution: 1bpm
- Measurement accuracy: +2bpm or 2% (The larger)

## Mechanical

- Alarm: SpO2 PR, Probe off, Finger out, Battery-low indicator
- Modes: Audio & Visual alarm, Information

## Data Management

- Transmission method: USB Data Cable
- ID setup: 1 – 127
- Data record: 72 hours

## Power Supply

- Type: (3) AA alkaline batteries or adapter
- Battery operation time: Work for 16 hours continuously at most
- AC adapter: Input Voltage – AC 100 – 240V, Input Frequency: 47 – 63 Hz, • • Output Current – 2A Max

Weight: 165 grams = 6oz

## Accessories:

- (3) AA alkaline batteries
- (1) Instructional manual
- (1) Adult finger probe
- (1) Water resistant poly bag
- The key for screwing the screw
- USB cable

Weight: 165 grams = 6 oz